

## Abutment Screw Driver Type and Torque Settings

Manufacturer	Torque N-cm	Driver for Titanium Screw	Driver for Gold Screw
Astra 3.0	15	0.050" -1.25mm hex	
Astra 3.5-4.0	20	0.050" -1.25mm hex	
Astra 4.5-5.0	25	0.050" -1.25mm hex	
Biohorizons® External, Internal, Internal Tapered	30	0.050" -1.25mm hex	
Biomet 3i™ External Hex	35	0.048" hex	Square
Biomet 3i™ Certain®	20	0.048" hex	
BlueSkyBio One Stage	30	Star & 0.048" hex	
BlueSkyBio Trilobe	30	Unigrip & 0.048" hex	
Branemark (Nobel Biocare)	35	Unigrip	Square
Camlog™ & Conelog®	20	0.050" -1.25mm hex	
Dentsply Xive / Friident	24	0.9mm hex	
Dentsply Ankylos® C/X	15	1mm hex	
Hiossen / Osstem HG mini (3.5)	20	1.2mm hex	
Hiossen/Osstem HG standard 4.0, 4.5, 5.0	30	1.2mm hex	
Imtec® 3M Endure™	35	0.050" -1.25mm hex	
Implant Direct™	30	0.050" -1.25mm hex	
Keystone/Lifecore Genesis & Prima™	30	Square	
Keystone/Lifecore Renova®	30	0.048" hex	
Keystone/Lifecore Restore®	30	0.048" hex	Square
Mega'Gen EZ Plus	35	0.048" hex	
Mega'Gen Rescue	45	0.048" hex	
Mega'Gen ExFeel internal (3.5)	25	1.2mm hex	
Mega'Gen ExFeel internal (4.1, 4.8, 5.5)	35	1.2mm hex	
Mega'Gen ExFeel external (3.3)	25	1.2mm hex	
Mega'Gen ExFeel external (3.75, 1.0, 4.5, 5.0, 5.5)	35	1.2mm hex	
MIS® Biocom & Seven	30	0.050" -1.25mm hex	
Neoss	30	Neoss Unigrip	
NobelActive™ & Replace™	35	Unigrip	
NobelActive™ 3.0	15	Unigrip	
OCO Biomedical	20	0.050" -1.25mm hex	
Southern Implant External	32	0.048" hex	Square & Unigrip
Southern Implant Tri Nex	32	Unigrip	
Southern Implant Octa	32	Star/Torx	
Straumann® Bone & Tissue Level	35	Star/Torx	
Sybron Pitt Easy	25	1.7mm hex	
Thommen SPI® 3.5	15	Four lobe	
Thommen SPI® 4.0, 4.5, 5.0, 6.0	20	Four lobe	
Zimmer	30	0.050" -1.25mm hex	

**CAD/CAM titanium and zirconia abutments, such as Atlantis™, 3i Bellatek™ and NobelProcera™, utilize the same torque and driver settings as the original implant manufacturer**