



**GNATHODONTICS, LTD**

— *Your Full Service Dental Laboratory* —

## **Crown Under Partial Impression Technique w/out sending partial**

### **Step One:** Standard Impression

- 1) Take a standard impression of the preparation site with the technique of your choice (without partial in place).
- 2) Take opposing impression and bite registration.

### **Step Two:** Make an index for relation of partial to the crown preparation (do not worry about the margins w/ this index, only pay attention to capturing the clasp(s) and rest(s)).

- 1) Inject BluDMousse Super Fast (recommended to reduce chair time) around the prep site, clasp(s) and rest(s).
- 2) Once the BluDMousse material has set (30 seconds), remove the partial with the index attached and trim the BluDMousse material so the outer portion of the clasp(s) and rest(s) are exposed.
- 3) Trim away any material that does not allow the index to be removed.
- 4) Remove the index from the partial by pushing index in an apical direction.

### **Step Three:** Temporize the preparation and send the impressions, index, bite registration, and lab slip to Gnathodontics.

### **Send the patient home wearing their partial**